Skyline Tutorial Webinar #8

DDA to Targeted: Differential Statistics with Skyline

With Brendan MacLean (Principal Developer)

Agenda

- Welcome from the Skyline team!
- DDA to Targeted: Differential Statistics
- Introduction with Brendan MacLean
 - Workflow and data set overview
- Tutorial with Brendan MacLean
 - DDA data processing review
 - Hypothesis generation from DDA data
 - Initial data review and Chorus Cloud Extraction
- Audience Q&A submit questions to Google Form:
 https://skyline.gs.washington.edu/labkey/qa4skyline.url



Chromatography-based Quantification

Hypothesis testing (Verification)

- SRM
- ▶ MS1 chromatogram extraction
- Targeted MS/MS (PRM)
- Data independent acquisition (DIA/SWATH)

Acquisition	Targeted	Survey
More Selective	PRM	DIA
Less Selective	SRM	MS1

Got HYPOTHESIS??





Multiple Instrument Vendors














Discovery to Targeted with Skyline





Case Study: ABRF iPRG 2014

	Fake Accession	Name	Origin	Molecular Weight
A	P44015	Ovalbumin	Chicken Egg White	45KD
В	P55752	Myoglobin	Equine Heart	17KD
C	P44374	Phosphorylase b	Rabbit Muscle	97KD
D	P44983	Beta-Galactosidase	Escherichia Coli	116KD
Е	P44683	Bovine Serum Albumin	Bovine Serum	66KD
F	P55249	Carbonic Anhydrase	Bovine Erythrocytes	29KD



Sample Preparation

	A	В	C	D	E	F	(fmol)
Sample 1	65	55	15	2	11	10	+ 200 ng yeast digest
Sample 2	55	15	2	65	0.6	500	+ 200 ng yeast digest
Sample 3	15	2	65	55	10	П	+ 200 ng yeast digest



Group Comparisons

	A	В	C	D	E	F	(fold change)
Sample 1-2	0.85	0.27	0.13	32.5	0.055	50	+ 200 ng yeast digest
Sample 1-3	0.23	0.036	4.33	27.5	0.91	1.1	+ 200 ng yeast digest
Sample 2-3	0.27	0.13	32.5	0.85	16.7	0.022	+ 200 ng yeast digest



Group Comparison Maxima

	A	В	C	D	E	F	(abs log2 fold change)
Sample 1-2	0.2	1.9	2.9	5.0	4.2	5.6	+ 200 ng yeast digest
Sample 1-3	2.1	4.8	2.1	4.8	0.1	0.1	+ 200 ng yeast digest
Sample 2-3	1.9	2.9	5.0	0.2	4.1	5.5	+ 200 ng yeast digest
Maximum	2.1	4.8	5.0	5.0	4.2	5.6	



DDA Runs Searched


	Identified yeast proteins
sample 1-a	3016
sample 1-b	3073
sample 1-c	2905
sample 2-a	2916
sample 2-b	2984
sample 2-c	2907
sample 3-a	2883
sample 3-b	2972
sample 3-c	2913

Comet, OMSA, MSGF+ - iProphet



Keys to Success with MS1 in Skyline

- Use Import DDA Peptide Search wizard
- Make sure you have ID annotations
 - Diagnose with Spectral Library Explorer
 - http://tinyurl.com/Skyline-missing-ids
- Review RT alignment in alignment viewer
- Got HYPOTHESIS??
- ▶ Review and manually adjust <5% of peaks
- Do the tutorial





Tutorial Discovery with Skyline

- ► Import DDA Peptide Search
- Data analysis

Four Processing Workflows Compared

- ▶ Peptides adjusted p value 0.05 remove single hits
 - 6 proteins + I false discovery (easily discounted)
- ▶ Peptides adjusted p value 0.01 remove single hits
 - 5 proteins (missing A)
- Proteins adjusted p value 0.05
 - 5 proteins (missing A) + 3 false discoveries (discounted)
- Proteins adjusted p value 0.01
 - ▶ 4 proteins (missing A & B) + 3 false discoveries (discounted)



Discovery versus Validation

- Discovery
 - asking your data what changed
- Validation
 - asking if there is evidence you candidates changed

Got HYPOTHESIS??



Learn More

- Webinars #1 (DDA) and #6 (Processing) tutorials
- Webinar #9:TBD
 - Tuesday, July 14
- Workshop in Rio de Janiero, August 31-September 2
- Workshop in Puerto Vallarta, November
- Weeklong Course at IIT-Bombay
 - December 10-14
- Weeklong Course in San Francisco
 - January?



Questions?

Ask any questions you have on differential statistics at the following form:

http://tinyurl.com/QA4Skyline

▶ Take the post-webinar survey:

http://tinyurl.com/Survey4Webinar



This ends this Skyline Tutorial Webinar.

Please give us feedback on the webinar at the following survey: http://tinyurl.com/Survey4Webinar

A recording of today's meeting will be available shortly at the Skyline website.

We look forward to seeing you at a future Skyline Tutorial Webinar.

