

Calibrated Quantification

Christopher M. Shuford

Laboratory Corporation of America® Holdings

LabCorp
Laboratory Corporation of America

Integrated ONCOLOGY
LabCorp Specialty Testing Group

National GENETICS INSTITUTE
LabCorp Specialty Testing Group

ViroMed LABORATORIES
LabCorp Specialty Testing Group

LithoLink STONE
LabCorp Specialty Testing Group

Endocrine SCIENCES
LabCorp Specialty Testing Group

LabCorp CLINICAL TRIALS

Monogram BIOSCIENCES
LabCorp Specialty Testing Group

Colorado COAGULATION
LabCorp Specialty Testing Group

Integrated GENETICS
LabCorp Specialty Testing Group

Dianon PATHOLOGY
LabCorp Specialty Testing Group


MedTox LABORATORIES
LabCorp Specialty Testing Group

ONE CHOICE
FOR ALL YOUR TESTING NEEDS


Why Calibrate?

*reduce magnitude of analytical variance
(and increase accuracy of relative quantification)*

Apolipoprotein E




Internal Calibration




$$\frac{[\text{NAT}]}{[\text{SIL}]} = \frac{I_{\text{NAT}}}{I_{\text{SIL}}} \longrightarrow [\text{NAT}] = \frac{I_{\text{NAT}}}{I_{\text{SIL}}} \times [\text{SIL}]$$

Internal Calibration Hierarchy



- 1 Full-length Protein:** Desiderio, D. M. and co-workers, *Biol. Mass Spectrom.* **1991**, 2(2), 149-156
Brun, V. and co-workers, *Mol. Cell. Proteomics* **2007**, 6(12), 2139-2149
- 2 Peptide Concatemer:** Beynon, R. J. and co-workers, *Nat. Methods* **2005**, 2(8), 587-589
Beynon, R. J. and co-workers, *Nat. Protocols* **2006**, 1(2), 1029-1043
- 3 (Partial) Peptide:** Barr, J. R. and co-workers, *Clin. Chem.* **1996**, 42(10), 1676-1682
Gygi, S. P. and co-workers, *P. Natl. Acad. Sci. USA* **2003**, 100(12), 6940-6945.

Internal Calibration




Pitfalls:


- Assume digestion is stoichiometrically complete (peptide and concatamers)
 - Results vary by digestion/denaturation conditions
- Reproducibility between digestions/days is critical
 - Internal Calibrator must be stable
- Response assumed to be linear

External Calibration: Multipoint


Qualified protein calibrator added to blank matrix



Without Internal Standard




With Internal Standard



External Calibration: Multipoint


Qualified protein calibrator added to blank matrix




Pitfalls:

- “Blank Matrix” must be essentially identical to sample matrix (requires verification)
- Qualified protein should have same structure as endogenous protein (size, conformation, PTMs, etc)
- External calibrators should be processed in parallel with unknown samples (added cost/time).
- External calibrators must be stable


External Calibration: Single-point




Without Internal Standard



With Internal Standard



External Calibration: Single-point



Pitfalls:

- “Blank Matrix” must be essentially identical to sample matrix (requires verification)
- Qualified protein should have same structure as endogenous protein (size, conformation, PTMs, etc)
- External calibrators should be processed in parallel with unknown samples (added cost/time).
- External calibrators must be stable
- ***Response assumed to be linear***

Additional Reading

Source of Analytical Variance in Proteomics

H.D. Cox *et al.*, *Clin. Chem.* **2016**, 60 (3), 541-548.

Peptide Degradation during Digestion

C.M. Shuford *et al.*, *Mol. Cel. Proteomics*, **2012**, 11 (9), 814-823.

More Peptides = More Confidence

Q. Fu *et al.*, *Clin. Chem.* **2016**, 62 (1), 198-207.

Single-point Calibration with Sample Pool

S.A. Agger *et al.*, *Clin. Chem.* **2010**, 56 (12), 1804-1813.

Analytical Verification in Translational Research

R.P. Grant and A.N. Hoofnagle, *Clin. Chem.* **2014**, 60 (7), 941-944.